

Cell Preparation Guide for Single Cell Sequencing Projects

Sample Quality

- ✓ minimal cell debris
- ✓ no cell aggregates
- ✓ at least 70% cell viability
- ✗ *dying cells increase ambient RNA contamination & cell clumping*

Sample Handling

- ✓ place samples on ice after resuspension/sorting ❄️
- ✓ store sorted/extracted samples not longer than **30 minutes** on ice before handing over the single cell core facility
- ✗ *prolonged sample handling negatively impacts sample quality*


Pipetting

- ✓ Speed: pipette cell suspensions slowly and gently
- ✗ *fast pipetting causes physical damage to cells by shearing forces*

Centrifugation conditions


- ✓ adjust centrifugation conditions by to sample type
 - **nuclei:** 500 rcf, 5 - 10 min
 - **small cells:** 300 rcf, 5 min
 - **large cells:** 150 rcf, 3 min
- ✓ avoid excessive centrifugation
- ✓ swinging bucket centrifuge is preferred over fixed angle

Centrifuges Classified by Rotor Type



Washing and Resuspension

- ✓ optimize buffer conditions: 1X PBS (Ca and Mg free) + 0.04% w/v BSA
(up to 1% w/v BSA) is recommended for most general sample preparation
⚠ for BEAM-Labeling: use 1X PBS + 2% FBS
- ✓ optimize buffer volumes, number of washes and centrifuge conditions to reduce cell loss and debris
- ✓ always leave behind ~50 µl supernatant to preserve pellet after centrifugation



- ✗ *discarding the entire supernatant during washing and resuspension steps may cause pellet disruption and massive cell loss*
- ✗ *buffers should not contain >0.1 mM EDTA, >3 mM magnesium or surfactants as these interfere with reverse transcription and GEM generation*

Straining and Filtering

- ✓ microfluidic channels of 10X Chromium X are <100 µm wide: use pore size of strainer that is larger than the maximum cell diameter, but small enough to catch larger clumps
- ✓ filter at the last wash step and not the final cell suspension